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HoMedics  
Service Center Dept. 168  
3000 Pontiac Trail  
Commerce Township, MI  
48390

Correo electrónico:  
[cservice@homedics.com](mailto:cservice@homedics.com)

# Instruction Manual

## Body Fat Analyzer and Scale

### How the Body Fat Analyzer and Scale Works:

This instrument uses BIA (Bioelectrical Impedance Analysis). BIA is one of the most accurate methods of measuring body fat and total body water. It sends a harmless signal through the fat and muscle in your body. This method calculates your personal weight, body fat and body water simultaneously, thus giving you a more accurate reading of your overall health and fitness. The body fat scale is also equipped with an **Athlete Mode** for athletes whose body build is different from non-athletes.

### Why it is Important to Measure Your Body Fat:

Body fat is a vital measure of your overall health and fitness. Research has shown that healthy body fat levels can help reduce your risk of developing serious health related problems, such as high blood pressure, high cholesterol, heart disease and diabetes.

The optimal Body Fat % (BF%) of an individual varies according to age and gender. Use the tables below as a guide:

**BODY FAT % RANGES FOR MEN<sup>†</sup>**

Rating	Age				
	20-29	30-39	40-49	50-59	60+
Low	<13	<14	<16	<17	<18
Normal	14-20	15-21	17-23	18-24	19-25
Moderate	21-23	22-24	24-26	25-27	26-28
High	>23	>24	>26	>27	>28

**BODY FAT % RANGES FOR WOMEN<sup>†</sup>**

Rating	Age				
	20-29	30-39	40-49	50-59	60+
Low	<19	<20	<21	<22	<23
Normal	20-28	21-29	22-30	23-31	24-32
Moderate	29-31	30-32	31-33	32-33	33-35
High	>31	>32	>33	>34	>35

<sup>†</sup> Source: University of Illinois Department of Food Science and Human Nutrition. Body Fat Percentage Calculator. [www.ag.uiuc.edu/~food-lab/ai/bfc.html](http://www.ag.uiuc.edu/~food-lab/ai/bfc.html)

**NOTE:**

Ideal body fat content is not the same for all people. Age, sex, and heredity are variables in this measurement. Always consult your physician to determine what is most ideal for you.

**Why it is Important to Monitor Your Total Body Water (TBW%) Percentage?**

Water is an essential component of the body and is one of the body's overall health indicators. Water makes up approximately 50-70% of the body's weight and is present proportionally more in lean tissue as compared to fat tissue.

Monitoring the level of body water can be a useful tool for one's health maintenance. Similar to body fat estimation, the TBW% (Total Body Water %) function provided in this scale is based on BIA. The estimated TBW% may vary according to your hydration level. For better accuracy, avoid fluctuations in hydration level prior to taking a measurement. The accuracy of the scale in estimating TBW% will also decrease with individuals suffering from diseases or illnesses that tend to accumulate water in the body. Please remember that the scale represents only a good approximation of your total body water.

The optimal TBW% of an individual varies according to age and gender. The table below may be used as a guide:

	BF % Range	Optimal TBW% Range <sup>††</sup>
Men	4 to 14%	70 to 63%
	15 to 21%	63 to 57%
	22 to 24%	57 to 55%
	25% and over	55 to 37%
Women	4 to 20%	70 to 58%
	21 to 29%	58 to 52%
	30 to 32%	52 to 49%

<sup>††</sup> Source: Derived from Wang & Deurenberg: "Hydration of fat-free body mass". American Journal Clin Nutr 1999; 69 833-841.

**Caution**

BIA (Bioelectrical Impedance Analysis) method determines your body fat and body water percentage by sending a harmless signal through the body. **Do not use this product if you have a pacemaker, artificial limbs, contraceptive devices, metal plates/screws or other internal medical devices. When in doubt, contact your physician.**

This analyzer will give accurate body fat and body water % readings for a majority of people, but is not intended for use by the following groups:

- Children:** Anyone under the age of 10 years
- Pregnant Women**

**Note:** This instrument is a personal monitor and should be used in a consistent manner for the most accurate reading. While readings might vary between this analyzer and other methods of measurement (hydrodensitometry or hand held calipers), your **changes** in body fat percentage shown by this instrument will be reflected accurately.

**Battery Installation**


1. Uses one 9-Volt alkaline battery (not included).
2. Remove Battery Compartment Cover on the bottom of the scale and insert the battery into the 9-Volt battery snap.
3. The scale needs to be initialized before first use or after battery replacement: Press lightly on the scale platform. The display will show "0000", then the scale will automatically turn off. Your scale is now ready for use.
4. Replace the battery when "Lo" is displayed.
5. Remove the battery if the product is not used for a long time.

## Getting Started

1. Select (Lb/in), (Kg/cm) or (St/in) by the slide switch on the bottom of the scale.

### Programming Personal Data into Memory

Press "Mode" button to enter Body Fat/Body Water mode.

1. **Set Height** - The height display blinks. Press the (▲)/(▼) button to adjust the height. Press "Mode" button to confirm.

2. **Set Age** - The age display blinks. Press the (▲)/(▼) button to adjust the age. Press "Mode" button to confirm.


3. **Set Gender** - The gender display blinks. Press the (▲)/(▼) button to adjust the gender. Press "Mode" button to confirm.

4. **Select Normal or Athlete** - The normal or athlete display blinks. Press the (▲)/(▼) button to adjust normal or athlete. Press "Mode" button to confirm.

**NOTE:** Athlete setting is for those aged 18 or above. This step is by-passed for those aged 17 and below.

• A "beep" sound indicates a successful setting. An "R" icon will be displayed.

5. **Select Memory Location** - Press the (▲)/(▼) button to select desired memory location (1-4). The memory number blinks. Press "Mode" to store the settings into the desired memory location. A "beep" sounds when the setting is complete. An "M" icon will be shown on the display.

**NOTE:** If the selected memory number is already occupied, a "■" icon will display in the lower right corner.

**Your Personal Data is Now Stored into Memory.**

## Instructions For Use


**Taking Your Body Fat and Body Water Reading:**

**You must have bare feet to make this measurement.**


To get the most accurate and consistent reading, wipe your feet with a damp cloth. Leaving them slightly moist, ensure the Analyzer/Scale is on a firm flat surface. The best time for using this analyzer is in the morning, before breakfast.

**To Recall Settings from Memory:**


1. While the power is off, press either the (▲) or (▼) button to toggle between Memory locations. Once you have selected the correct memory location, Press "Mode" to confirm.
2. "R" icon will be displayed on the screen. The screen displays the height, age, gender and normal/athlete selection on the toggle screen sequentially after 2 seconds.


- Step on the scale with bare feet. Stand still while the unit measures your weight. The weight will display for 2 seconds.


- Then it starts to measure body fat and total body water. The measurement pattern will be shown on the display during these measurements.

- The measurements finish with 1 long "beep". The fat percentage, total body water and weight results will toggle in 2 second intervals. The unit will turn off automatically after 8 seconds.

The "Status Bar" at the bottom of the screen will show your general condition. See "General Information About Body Fat" for recommended body fat levels.

#### To Clear Memory:

- While the power is off, press either the (▲) or (▼) button to toggle between Memory locations.
- If the selected memory is occupied, press both (▲)/(▼) at the same time to clear the content of the selected Memory.


#### To Turn Off the Scale:

To turn the scale off, press and hold (MODE) key for 2 seconds. Or, if no key has been pressed, the unit turns off automatically after 15 seconds.


#### Body Fat & Body Water Measuring Operation without Memory:

- Follow steps 1-4 in "Programming Personal Data into Memory" section

- Step on the scale with bare feet.** The unit will now calculate your weight reading.

- Stay on the scale. Stand still while the unit measures your body fat and total body water %. The measurement pattern will be shown on the display during these measurements.


- The BFA measurements finish with 1 long "beep". The fat percentage, total body water and weight results will toggle in 2 seconds intervals. The unit will turn off automatically after 8 seconds.

#### Weight Ready Only:

- Step onto the scale; stand still while the scale computes your weight.

2. Your weight will be shown on the LCD display as shown below.

3. The scale turns off by itself after use.

**NOTE:** Place the Analyzer/Scale on a firm, flat surface to achieve best accuracy.

### TROUBLE SHOOTING

1. **You must have bare feet when using body fat and total body water measurements.**  
To get the most accurate and consistent reading, wipe your feet with a damp cloth, leaving them slightly moist before stepping on scale.
2. Be sure your feet are centered on the electrodes.
3. **Display Err0** Initialization error. The user is advised to take the measurement again for an accurate result.
4. **Display Err1** Wrong usage or body not stable. Repeat measurement again.
5. **Display Err2** Overload Warning. Remove the weight immediately; otherwise, permanent damage to the scale will occur.
6. **Display Err3** Contact error. Impedance cannot be measured. Make sure that you are standing still on the scale and maintain maximum contact between your feet and metal contacts. Clean bottom of feet with a damp cloth, leave slightly moist, and repeat measurement again.
7. **Display Err4** Body fat percentage is beyond the technical limit. Clean bottom of feet with a damp cloth, leave slightly moist, and repeat measurement again.
8. **Display Err5** Percentage total body water is out of range. Moistening your feet may help to improve the contact.
9. **Display Lo** Replace the battery when "Lo" is displayed.
10. The condition of the skin on the bottom of your feet can affect the reading. The natural effects of aging or activity can make this skin hard. Take the reading with clean, slightly moist feet for best accuracy.

### Care and Maintenance

1. Do not disassemble the scale other than replacing the battery; it contains no user serviceable parts. Damage to the scale may occur as a result of improper handling.
2. Remove battery when the scale is not used for a prolonged period of time.
3. Clean the scale after use with a lightly dampened cloth. Do not use solvents or immerse the unit in water.
4. Avoid excessive impact or vibration to the scale, such as dropping it onto the floor.


The product is intended for home/consumer use only; it is not intended for professional use in hospitals or medical facilities.

**Note: If you are having any problems operating this Body Fat Analyzer, please call 800-466-3342, between the hours of 8:30 a.m. to 5:00 p.m. EST**

### Specifications

- Uses BIA (Bioelectrical Impedance Analysis) to determine body fat and total body water percentage
- High precision STRAIN GAUGE technology to determine weight
- Weight Capacity: 330 Lb, 150 Kg or 23 St
- Weight Graduation: 0.2 lb, 100 g, 1/4 lb (Dashes shown below weight reading in Stone mode)
- Body Fat Graduation: 1.0% (in the range from 4% - 50%)
- Body water graduation: 1.0%
- Height Range: 3'3.4" - 7' 2 " (100.0 - 220 cm)
- Age Range: 10 - 99 years
- Measures weight, body fat and total body water simultaneously in about 10 seconds
- Athlete or Normal mode selection
- Weight-only feature
- Auto on/auto off function for convenience and conserving power
- Data can be stored for up to 4 individuals
- User selectable units of measure (Lb/in), (Kg/cm), or (St/in)
- Low-battery indicator
- Operates on one 9-Volt alkaline battery (not included)

### LCD Display

## LIFETIME WARRANTY

HoMedics, Inc., guarantees this product free from defects in material and workmanship for the life of the original purchaser from the date of original purchase, except as noted below.

This HoMedics product warranty does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorized accessory; alteration to the product; or any other conditions whatsoever that are beyond the control of HoMedics. This warranty is effective only if the product is purchased and operated in the USA. A product that requires modification or adaptation to enable it to operate in any country other than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under warranty. HoMedics shall not be responsible for any type of incidental, consequential or special damages.

To obtain warranty service on your HoMedics product, either hand deliver or mail the unit and your dated sales receipt (as proof of purchase), postpaid, along with check or money order in the amount of \$5.00 payable to HoMedics, Inc. to cover handling.

Upon receipt, HoMedics will repair or replace, as appropriate, your product and return it to you, postpaid. If it is appropriate to replace your product, HoMedics will replace the product with the same product or a comparable product at HoMedics' option. Warranty is solely through HoMedics Service Center. Service of this product by anyone other than HoMedics Service Center voids warranty.

This warranty provides you with specific legal rights. You may have additional rights which may vary from state to state. Because of individual state regulations, some of the above limitations and exclusions may not apply to you.

**For more information regarding our product line in the USA, please visit: [www.homedics.com](http://www.homedics.com)**

**For service or repair - do not return this unit to the retailer. Contact HoMedics Consumer Relations at 1-800-HOMEDIC (1-800-466-3342).**

**Mail To:**

**HoMedics  
Service Center Dept. 168  
3000 Pontiac Trail  
Commerce Township, MI  
48390**

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[www.ag.uiuc.edu/~food-lab/ai/bfc.html](http://www.ag.uiuc.edu/~food-lab/ai/bfc.html)

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