BW Series Compact Bench Scale

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Balance de table compacte, série BW

BW SERIES COMPACT BENCH SCALE EN-1

| | |
|---|-------|
| 1. INTRODUCTION | EN-2 |
| 1.1 Product Description..... | EN-2 |
| 1.2 General Features..... | EN-2 |
| 1.3 Safety Precautions | EN-2 |
| 2. INSTALLATION | EN-3 |
| 2.1 Unpacking | EN-3 |
| 2.2 Installing Components..... | EN-3 |
| 2.3 Selecting the Location | EN-3 |
| 2.4 Leveling the Scale..... | EN-3 |
| 2.5 Connecting Power..... | EN-4 |
| 2.5.1 AC Power..... | EN-4 |
| 2.5.2 Battery Power..... | EN-4 |
| 2.6 Initial Calibration | EN-4 |
| 3. OVERVIEW OF CONTROLS AND FUNCTIONS | EN-5 |
| 4. OPERATIONS..... | EN-6 |
| 4.1 Switching the Unit On and Off..... | EN-6 |
| 4.2 Manual Tare..... | EN-6 |
| 4.3 Zero Operation..... | EN-6 |
| 4.4 Changing Weighing Units | EN-6 |
| 5. SCALE SETTINGS | EN-7 |
| 5.1 Scale Increment..... | EN-7 |
| 5.2 Backlight | EN-7 |
| 6. CALIBRATION AND SEALING | EN-8 |
| 6.1 Calibration | EN-8 |
| 6.2 Legal for Trade Sealing..... | EN-10 |
| 7. MAINTENANCE | EN-11 |
| 7.1 Cleaning..... | EN-11 |
| 7.2 Troubleshooting..... | EN-12 |
| 8. TECHNICAL DATA..... | EN-13 |
| 8.1 Specifications..... | EN-13 |
| 8.2 Drawings..... | EN-14 |
| 8.3 Compliance | EN-15 |
| WARRANTY | EN-16 |

1. INTRODUCTION

This manual contains installation, operation and maintenance instructions for the BW Series Compact Bench Scale. Please read the manual completely before installation and operation.

1.1 Product Description

The BW Series Compact Bench Scale is an economical portioning, packing and general weighing scale designed to provide fast and accurate weighing performance in washdown environments. With an internal sealed rechargeable lead-acid battery, IP65 housing and high-speed display update, the BW Series Compact Bench Scale offers portability, durability and productivity in one practical design.

1.2 General Features

- Capacity models: 1.5kg/3lb, 3kg/6lb, 6kg/15lb, 15kg/30lb
- Single and dual display models
- Default and LFT displayed resolution 1:3,000
- High-profile ABS housing, stainless steel pan (pan is dishwasher-safe)
- IP65 dust and water jet protection
- 4-adjustable leveling feet with integrated front-panel leveling indicator
- Bright red LED weight display: 6-digit, 7-segment, 0.55"/14 mm high characters
- 4-button control panel: On/Units Off, Check, Tare, Zero
- 3-weighing units: g, oz, lb
- Battery charge indication
- AC power or 15-hour rechargeable lead-acid battery
- Over-Under checkweighing function with audible alarm

1.3 Safety Precautions

For safe and dependable operation of this scale, please comply with the following safety precautions:

- Verify that the input voltage printed on the AC Adapter and the plug type matches the local AC power supply.
- Make sure that the power cord does not pose an obstacle or tripping hazard.
- Disconnect the scale from the power supply when cleaning the scale.
- Do not operate the scale in hazardous or unstable environments.
- Do not drop loads on the platform.
- Use only approved accessories and peripherals, as available.
- Operate the scale only under ambient conditions specified in these instructions.
- Service should be performed by authorized personnel only.

2. INSTALLATION

2.1 Unpacking

Unpack and verify that the following components have been included:

- BW Scale Unit
- Weighing Pan
- AC-DC Adapter
- Instruction Manual


Save the packaging material. This packaging ensures the ideal protection for the storage or transport of the product.

2.2 Transport Protection

BEFORE USING THE SCALE:

Remove the 3 shipping protection screws on the spider sub-platform, then remove the 3 plastic column inserts under the spider. Re-install the protective components if shipping the scale to ensure proper transport protection.


2.3 Selecting the Location

Operate the scale on a firm, level surface. Avoid locations with rapid temperature changes, excessive air currents, vibrations, electromagnetic fields, heat or direct sunlight.

2.4 Leveling the Scale

Adjust the leveling feet* until the bubble is centered in the circle of the level indicator (located on the front panel). Ensure that the scale is level each time its location is changed.

* **IMPORTANT:** To allow proper air flow and minimize condensation inside the scale, unscrew the 2 rear leveling feet slightly, then level the scale accordingly for normal operation. However, when cleaning the scale, all 4 leveling feet must be screwed in tightly to avoid entry of water inside the scale.

2.5 Connecting Power


2.5.1 AC Power

Verify that the intended AC power source matches the AC adapter rating. Unplug* the rubber access cap underneath the scale, connect the adapter to the power receptacle through the cut-out, then plug the AC adapter into a power outlet.

* **NOTE:** When cleaning the scale, turn off the unit, remove the AC adapter and plug the rubber access cap securely to avoid entry of water inside the scale.

2.5.2 Battery Power

The scale will begin charging with the AC adapter connected accordingly, and automatically switches to battery operation if there is a loss of power. Normal operating time on a full charge is 50 hours on a 15-hour charge cycle.

Low battery charge is indicated by the flashing LED over the  indicator. AC power or a fully charged battery is indicated by the illuminated LED over the  indicator. A blinking  indicates the battery is charging.

NOTES:

- DO NOT continuously charge the scale longer than 24 hours.
- To maximize battery power, the scale can be set up to enter Sleep mode (only decimal points are displayed) after 2 minutes of inactivity, or automatically turn off after 15 minutes of inactivity (Section 5.2).
- Charging the scale must be performed in a dry environment.
- Plug the rubber access cap to the adapter receptacle to ensure the scale is protected when operated in a wet environment.
- If the scale is not used within 3 months, recharge and discharge the battery fully 3 times to maximize the battery condition.
- Dispose of the lead-acid battery according to local laws and regulations.

CAUTION


BATTERY IS TO BE REPLACED ONLY BY AN AUTHORIZED OHAUS SERVICE DEALER. RISK OF EXPLOSION CAN OCCUR IF REPLACED WITH THE WRONG TYPE OR CONNECTED IMPROPERLY.

2.6 Initial Calibration

When the scale is operated for the first time, a Span Calibration is recommended to ensure accurate weighing results. Before performing the calibration, be sure to have the appropriate calibration weights. Refer to Section 6 for calibration procedures.

3. OVERVIEW OF CONTROLS AND FUNCTIONS


- 
 - Turns the scale on
 - Long-press turns the scale off
 - Short-presses display successive units of measure
 - In SETUP Mode, “Yes” accepts/confirms the displayed setting
 - In CHECK Mode, accepts/confirms the displayed value

- 
 - Inputs the weight of the object on the weighing pan as a Tare value
 - In SETUP Mode, “No” toggles through the available settings
 - In CHECK Mode, increments the value of the selected digit

- 
 - Zeros the display
 - With the scale off, long-press together with **ON/UNITS OFF** enters User Setup Mode
 - In SETUP Mode, “Exit” ends setup and exits to “Save”
 - In CHECK Mode, ends Over/Under setup and exits to weighing mode

- 
 - Initiates CHECK Over/Under setup
 - In SETUP Mode, “Back” reverts to the previous parameter
 - In CHECK Mode, shifts to the next digit

- 
 - LED annunciator indicates Center of Zero

- 
 - LED annunciator indicates battery is fully charged (blinking if charging)

- 
 - LED annunciator indicates displayed weight is less than the preset Under value

- 
 - LED annunciator indicates displayed weight is more than the preset Over value

- 
 - LED annunciator blinks to indicate the battery voltage is low

- 
 - LED annunciator indicates Net (Gross minus Tare) value

- 
 - LED annunciator indicates the active weighing unit

4. OPERATIONS

4.1 Turning the Unit On and Off

With the scale off, press **ON/UNITS OFF** to turn on. Make sure the weighing pan is empty before turning the scale on. Allow 15-30 minutes for the scale to warm up before use. With the scale on, long-press **ON/UNITS OFF** to turn off.

4.2 Manual Tare

Place a container on the weighing pan, then press **TARE**. The display will show a net value of "0" as the Tare weight is registered. Remove the container, and the resulting negative value displayed will be equivalent to the Tare weight. To clear the Tare value, press **TARE** with the pan empty.

4.3 Zero Operation

Press **ZERO** to zero the stable weight display (the weight must be within $\pm 2\%$ of scale capacity). Center of Zero is indicated by the LED annunciator.

4.4 Checkweigh Target Setting

Compares the weight of an item against a pre-set Over-Under range.

To activate Checkweighing Mode, refer to user setup (Section 5.2).

- Press **CHECK** to initiate entering Over-Under parameters.
- The value of the "OVER" target limit will be displayed first, with the active digit flashing (the "OVER" LED annunciator is lit).
- Enter the desired "OVER" limit value by pressing **-No** to increment the active digit or **-Back** to shift to the next digit.
- Press **-Yes** to accept the displayed value and proceed to the "UNDER" target limit (the "UNDER" LED annunciator is lit).
- Enter the desired "UNDER" limit value by pressing **-No** to increment the active digit or **-Back** to shift to the next digit.
- Press **-Exit** to accept the displayed value and exit into checkweighing mode.

NOTES:

- "OVER" or "UNDER" LED annunciators will light if the weight of an item placed on the pan is over or under the pre-set limits, respectively. "OVER" and "UNDER" LED annunciators will light simultaneously if the weight is within the pre-set limits.

- The LED indication is complemented by an audible alert feature (Sec. 5.2).
- Over-Under values are retained when the unit is turned off.
- Over-Under values are not retained when the weighing unit or resolution setup are changed.

5. SCALE SETTINGS

5.1 Accessing the User Setup Mode

With the scale on, press **ON/UNITS OFF** and **CHECK** simultaneously. **"SETUP"** then **"S1 OFF"** is displayed. To navigate in the User Setup Mode:

- Press **ON/UNITS OFF-Yes** to accept the displayed setting
- Press **TARE-No** to toggle through available settings
- Press **ZERO-Exit** to end setup and proceed to "Save"
- Press **CHECK-Back** to revert to the previous parameter

5.2 User Setup Parameters

| Parameters | Description | Settings |
|-------------|---|---|
| S1 | Scale automatically turns off after 15 minutes of inactivity. | ON - enables this function. OFF -disables this function. |
| S2 | Scale enters sleep mode after 2 minutes of inactivity. | ON - enables this function. OFF -disables this function. |
| S3 | CHECK activation and Alert Range | 0 - Check mode is off. <i>Scale beeps when displayed weight is:</i> 1-...within the preset Over and Under values 2-...outside the preset Over and Under range 3-...within the preset Over and Under values (Over-Under LED's do not light within preset range) 4-...outside the preset Over and Under range (Over-Under LED's do not light within preset range) 5- Alarm is off (Over-Under LED's do not light within preset range) 6- Alarm is off |
| S4 | Averaging Level | 0, 1, 2 , 3 (most filtering, longer stabilization time) 4, 5, 6, 7 (corresponds to 0, 1, 2, 3 but with slower display update) |
| S5 | Power-on Weighing Unit | lb , oz, g |
| S6 | LED Brightness Level | 0 ,1,2,3 (dimmmst). |
| dsp | Display Type | OFF - displays weight continuously. ON - displays first stable weight. |
| Save | Save or Abort | Save - saves changes and reverts to weighing mode. Abort - reverts to weighing mode without saving changes. |

Note: Bold settings are defaults.


6. CALIBRATION AND SEALING

6.1 Calibration

For best results, calibrate the scale at regular intervals. Temperature changes, geographic gravity variations, altitude changes and abuse are few reasons why a scale may need recalibration.

Note: The following procedure requires the use of hand tools and calibration weights. If unable to perform the calibration, contact an authorized Ohaus service organization.

1. Remove the battery cover to access the two calibration pins located in a slot at the bottom of the scale. With the scale on under battery power, short the calibration pins (a small slotted screwdriver may be used).

Note: When accessing the bottom of the scale, avoid placing the scale upside down on the pan or sub-platform. Place the scale on its side.

2. **SETUP** then the first parameter **CIY** is displayed.
3. Press **-Yes** successively to move through the succeeding parameters: **dEF NO**, **S1**, **S2**, **S3**, **S4**, **S5**, **S6**, **dSP**, **S8**, **S9**, **S10**, **S11**, **S12**, **bAt**, **GEO**.

Note: Do not change the settings of these parameters. These are service setup parameters intended for authorized servicing use only.

4. After **GEO**, **CAL NO** is displayed. Press **-No** to switch the setting to **CAL YES**.
5. **CAL YES** is displayed. Press **-Yes** to enter calibration sequence.

BW SERIES COMPACT BENCH SCALE EN-9

6. "g" is displayed. This is the Calibration Weight Unit. Press **-Yes** to accept (or press **-No** to switch to "lb", then press **-Yes**) and move on to the next parameter. Refer to the following chart for the next steps.

Selected calibration weight unit: g

| | | | | |
|-------------------|---------------|-------------|--------------|---------------|
| Scale Capacity | 1.5 kg / 3 lb | 3 kg / 6 lb | 6 kg / 15 lb | 15 kg / 30 lb |
| Displayed Setting | 1500 | 3000 | 6000 | 15000 |
| 2/3 of Capacity | 1000 / 1 kg | 2000 / 2 kg | 4000 / 4 kg | 10000 / 10 kg |
| Full Capacity | 1500 g | 3000 g | 6000 g | 15000 g |

Selected Calibration weight unit: lb

| | | | | |
|-------------------|---------------|-------------|---------------|---------------|
| Scale Capacity | 3 lb / 1.5 kg | 6 lb / 3 kg | 15 lb / 6 kg | 30 lb / 15 kg |
| Displayed Setting | 3000 | 6000 | 15000 | 30000 |
| 2/3 of Capacity | 2000 / 2 lb | 4000 / 4 lb | 10000 / 10 lb | 20000 / 20 lb |
| Full Capacity | 3000 / 3 lb | 6000 / 6 lb | 15000 / 15 lb | 30000 / 30 lb |

7. "xxxx" is displayed. This is the Full Scale Capacity of the scale. Make sure the displayed value matches the scale model capacity then press **-Yes** to accept and move on to the next parameter.
8. "-----" is displayed. Make sure the pan is empty and stable at this time then press **-Yes**. The display will count down from 5 to 0 as the scale zero is registered.
9. "xxxx" is displayed. This is the Calibration Point corresponding to 2/3 of the Full Scale Capacity value. Place the indicated calibration weight on the center of the pan. With the scale stable, press **-Yes**. The display will count down from 5 to 0 as the scale registers the calibration weight.
10. "xxxx" is displayed. This is the Calibration Point corresponding to the Full Scale Capacity value. Place the indicated calibration weight on the center of the pan. With the scale stable, press **-Yes**. The display will count down from 5 to 0 as the scale registers the calibration weight.
11. "SAVE" is displayed when calibration was registered. Remove the calibration weight, then press **-Yes**. The scale will display a self-test then return to normal weighing mode. Or, press **-No** to abort calibration without saving any changes.


Note: At any time during the calibration process, long-press **ON/UNITS OFF** to turn off the scale and abort the process without saving any changes.

6.2 Legal for Trade Sealing

After the scale has been found to comply with local regulations by relevant approving authorities, the scale may be sealed in the following manner:


Lead and Wire Seal.

Secure sealing plate with through-hole screws. Run the wire seal into the through-holes then press the lead seal in place.

Sticker Seal.

Secure sealing plate with regular pan-head screws then apply the sticker over one of the screws.

7. MAINTENANCE

7.1 Cleaning

Comply with safety precautions when cleaning and maintaining the scale

- Turn off and disconnect the scale from the power supply during cleaning.
- **DO NOT** open the scale; no user repairable components inside. Service should be performed by authorized personnel only.
- **DO NOT** place the scale upside down. If accessing the bottom, place the scale on its side.
- **DO NOT** carry the scale by the pan or sub-platform. Use the handholds on the sides of the scale housing.
- **DO NOT** allow liquids to enter the scale.
- **DO NOT** use high pressure when cleaning the scale.
- Inspect the scale regularly. Ensure that the housing, pan and sub-platform are clean and free from foreign material.
- Wipe the exterior components of the scale using a soft cloth dampened with water and a mild detergent (**DO NOT** use acids, alkalis, strong solvents or abrasive cleaning agents).
- Screw in all 4 leveling feet tightly to avoid entry of water inside the scale during cleaning. However, when the scale is put back into operation, remember to loosen the two rear feet (and adjust leveling accordingly) to allow proper air flow and minimize condensation inside the scale.
- Plug the rubber cap to the adapter receptacle underneath the scale to avoid entry of water inside the scale.
- The pan is dishwasher safe and may be placed with normal loads. All other parts must not be placed in a dishwasher.

7.2 Troubleshooting

The following table lists common problems, possible causes and remedies. If the problem persists, contact Ohaus or your authorized Ohaus dealer.

| Symptom | Possible Causes | Remedy |
|------------------------------------|---|---|
| Scale will not turn on | AC power not connected to scale or to outlet. No power in outlet. Battery discharged. | Connect scale to power. Connect scale to power and charge battery. Connect to different power outlet. Replace battery. |
| Battery fails to charge fully. | Battery is defective or past its useful life. | Have the battery replaced by an authorized Ohaus service dealer. |
| Weight reading does not stabilize. | Unstable environment. Interference under pan or sub-platform. | Ensure a stable environment. Ensure that the weighing pan and sub-platform are unobstructed. |
| Scale does not display accurately. | Improper calibration | Calibrate the scale using proper calibration weights. |
| "nnnnn" displayed | Overload condition: Load on the pan exceeds scale capacity +9d | Reduce the load on the pan. Re-zero the scale. |
| "uuuuu" displayed | Underload condition: No load on the pan, exceeds scale zero -9d | Install the pan on the sub-platform. Re-zero the scale. |
| "- - - - -" displayed | Unstable condition. Cannot establish zero. | Situate in a stable environment. Relieve internal pressure from the scale by loosening the leveling feet. Clear possible obstructions on the vent holes under the rear feet. Remove the weight on the pan that may exceed the power-on zero range. |
| "E 11" displayed | RAM error | Turn the scale off and on again. |
| "E 16" displayed | ROM error | Recalibrate the scale. |
| "E 18" displayed | EEPROM error | Return scale for service |
| "E 48" displayed | Checkweigh target setting error | Review setup and re-enter correct Over-Under values |
| "Error" displayed | Software running error | Turn the scale off and on again. Replace battery. |



8. TECHNICAL DATA

8.1 Specifications

| Model * | BW1.5US | BW3US | BW6US | BW15US |
|------------------------------------|--|---|---|---|
| Capacity x Readability | 1.5kg x 0.5g 3lb x 0.001lb 50oz x 0.02oz | 3kg x 1g 6lb x 0.002lb 100oz x 0.05oz | 6kg x 2g 15lb x 0.005lb 200oz x 0.1oz | 15kg x 5g 30lb x 0.01lb 500oz x 0.1oz |
| Displayed / LFT Resolution (lb/kg) | 1:3,000 | | | |
| Construction | ABS plastic housing with stainless steel pan | | | |
| Protection | IP65 | | | |
| Pan Dimensions | 9"W x 7"D / 230W x 180D mm | | | |
| Leveling Components | 4 x adjustable feet with integral panel mounted level bubble indicator | | | |
| Display * | 1 x red LED; 6-digit, 7-segment, 0.55"/14mmH characters | | | |
| Weighing Units | lb, g, oz | | | |
| Keyboard | 4 x tactile mechanical keys: ON/UNITS OFF, CHECK, TARE, ZERO | | | |
| Calibration Weight | 2/3 and Full Scale Capacity, lb/g | | | |
| Over Range Capacity | Full Scale Capacity +9d | | | |
| Tare Range | Full Capacity by Subtraction | | | |
| Zero Range | Power-on: $\pm 10\%$ / Keyboard: $\pm 2\%$ Full Scale Capacity | | | |
| Stabilization Time | within 2 seconds | | | |
| Power | 12VDC 0.84A Max AC Adapter / 6V 5Ah rechargeable lead acid battery | | | |
| Battery Life | 50 hours continuous use with 15 hour recharge time | | | |
| Safe Overload | 150% Full Scale Capacity | | | |
| Operating Temperature Range | 10°F to 104°F / -10°C to 40°C | | | |
| Storage Temperature Range | 0°F to 122°F / -25°C to 50°C | | | |
| Humidity Range | 10% to 85% humidity, non-condensing | | | |
| Scale Dimensions | 9"W x 12.2"D x 5"H / 230W x 310D x 127H mm | | | |
| Shipping Dimensions | 11.6"W x 15"D x 7.5"H / 295W x 381D x 190H mm | | | |
| Scale Unit Weight | 9.7 lb / 4.4 kg | | | |

* Note: BWxxTUS = dual display model

8.2 Drawings

8.3 COMPLIANCE

Compliance to the following standards is indicated by the corresponding mark on the product.

| Marking | Standard |
|---|---|
|  | This product conforms to the EMC directive 2004/108/EC, the Low Voltage Directive 2006/95/EC and the Non-automatic Weighing Instruments Directive 90/384/EEC. The complete Declaration of Conformity is available from Ohaus Corporation. |
|  | UL60950-1: 2003 |

EC Emissions Note

This device complies with EN55011/CISPR 11 Class B Group 1.

| | |
|--|--|
|   | <p>Disposal</p> <p>In conformance with the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE) this device may not be disposed of in domestic waste. This also applies to countries outside the EU, per their specific requirements.</p> <p>Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment.</p> <p>If you have any questions, please contact the responsible authority or the distributor from which you purchased this device.</p> <p>Should this device be passed on to other parties (for private or professional use), the content of this regulation must also be related.</p> <p>Thank you for your contribution to environmental protection.</p> <p>For disposal instructions in Europe, refer to www.ohaus.com/weee.</p> |
|--|--|

FCC Note

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

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ISO 9001 Registration

In 1994, Ohaus Corporation, USA, was awarded a certificate of registration to ISO 9001 by Bureau Veritas Quality International (BVQI), confirming that the Ohaus quality management system is compliant with the ISO 9001 standard's requirements. On May 15, 2003, Ohaus Corporation, USA, was re-registered to the ISO 9001:2000 standard.

LIMITED WARRANTY

Ohaus products are warranted against defects in materials and workmanship from the date of delivery through the duration of the warranty period. During the warranty period, Ohaus will repair, or, at its option, replace any component(s) that proves to be defective at no charge, provided that the product is returned, freight prepaid, to Ohaus.

This warranty does not apply if the product has been damaged by accident or misuse, exposed to radioactive or corrosive materials, has foreign material penetrating to the inside of the product, or as a result of service or modification by other than Ohaus. In lieu of a properly returned warranty registration card, the warranty period shall begin on the date of shipment to the authorized dealer. No other express or implied warranty is given by Ohaus Corporation. Ohaus Corporation shall not be liable for any consequential damages.

As warranty legislation differs from state to state and country to country, please contact Ohaus or your local Ohaus dealer for further details.

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P/N 80252763 A © Ohaus Corporation 2008, all rights reserved / todos los derechos reservados / tous droits réservés