

How To Use

4 STEPS TO GET IN TUB



Approach **step2tub**, grab firmly onto each bar with both hands. Please be sure to balance yourself prior to stepping up on to the base.



With both hands firmly on the grab bars gently step up onto the base.



While still balancing on both grab bars evenly, turn to the side and prepare to step into the tub.



Lift one leg first over the side base of the tub and step gently into the tub. Repeat this step for your other leg and fully step into the bathtub.

5 STEPS TO GET OUT OF TUB



After you have thoroughly dried off and the area around the **step2tub** is free of water, place your hands firmly on both grip bars.



With both hands firmly on the grip bars, lift each leg one at a time to step out of the tub.



Come to a standing-still position with both feet firmly on the base of the step.




(When ready) Turn slightly to the side with both feet parallel to one another. Be sure to continue to place both hands firmly on either side of the grab bars.



With both hands still gripped firmly on either side of the grab bars, gently step off the base and on to the floor. Please ensure that the area around the **step2tub** is free of any water.

1. Place all parts on a flat working surface. Remove all plastic wrap, corner step protectors and foam wrapping before assembly.



2. Slide the legs into the base of the step stool making sure the brass tension tabs lock into place.



3. Repeat the process with all the four legs and set the step stool upright.



4. Prior to inserting the two handles into the legs, remove the circular black end pieces (as shown in picture below) on just the back side of the grab bars where the rounded clamp bars will slide into (you only need to remove a total of two of those pieces). Then, slide the rounded clamp bars (facing downwards towards the floor) onto each grab bar in an upward motion. Make sure these round-shaped bars are facing the back of the unit. Next, place the circular black end pieces (the two you just removed) back onto the bottom of the grab bars just below the rounded-clamped bars.

Now, insert the two handles into the legs making sure the handles are the same height by pushing in the brass tab adjustments. Please ensure that the holes are aligned properly on both sides as the holes of each bar need to be facing outwards.


NOTE:

Make sure they are locked in place and aligned properly before using.



5. To adjust height, press in each brass tab on each leg to adjust the height of the leg making sure all four legs are the same height.

The legs are adjustable to three positions by 1" increments.



WARNING:

Never use the step2tub unless the brass tabs are securely locked in place to avoid serious injury. Always check the brass tabs prior to using the step stool.

6. Customize height of the step2tub in increments of 1" and the handle by 1/2".

*Please note that the holes on both grab bars need to be facing outwards for proper use.



7. Apply both bathtub grab bars to the bathtub using the step2tub to properly measure and determine the correct alignment. Once they are aligned, tighten both bars as firmly as possible to ensure maximum stability. Twist handle to tighten by hand.



8. Lower the rounded bars on top of the bathtub grab bars and push the poles firmly into place. Please be sure to double check that the unit is firmly attached to the bathtub. Your safety is our priority!

