

I N S T R U C T I O N S



No-Wire™ Walker Basket (Assembled)

Catalog Number 70317-0001


U.S. Patent D451,054

Other Patents Pending

Installation Instructions (Can be installed on the inside or outside of the walker.):

1. Hang the walker basket on the cross bar at the front of the walker as shown below. If there is more than one cross bar, hang the basket on the lowest one.
2. Slide the retractable arms out from under the basket and snap on to the legs of the walker to keep the basket stable and upright.

NOTE: On walkers where there are features that prevent the basket from hanging properly, hook and loop strips can be used to secure the basket to the cross bar. Hook and loop can be purchased at your local craft, home healthcare or hardware store.



Maddak Inc.

6 INDUSTRIAL ROAD • PEQUANNOCK, NJ 07440-1993 USA

973-628-7600 • fax 973-305-0841 • email custservice@maddak.com • www.maddak.com