

FAT LOSS monitor



HBF-306C

*Accurate AND Consistent
Readings in Seconds!*

Reach Your Fitness Goals Faster!

- ✓ Displays body fat % and body mass index (BMI)
- ✓ Uses bioelectric impedance to measure body fat versus lean body weight
- ✓ Portable, easy to use

Features:

- ✓ Stores 9 personal profiles plus guest
- ✓ Measures in normal or athlete mode based on the Omron® FIT index
- ✓ Set to US or Metric system
- ✓ Instruction manual includes recommended measurement times, BMI and Body Fat % guidelines by age and gender, how to calculate your FIT index and more!
- ✓ One year limited warranty



FAT LOSS monitor



PACKAGING SPECIFICATIONS:

CARTON:	QTY:	APPROX. WEIGHT:	APPROX. DIMENSIONS: (w x h x d)
Unit:	1 Each	0.84 lbs.	8.66" x 2.3" x 6"
Inner:	3 Each	3.06 lbs.	7.8" x 6.5" x 9.3"
Master:	12 Each	14.39 lbs.	15.9" x 14.1" x 10.6"
UPC Code:	0 73796 30630 4		

UNIT SPECIFICATIONS:

Power Supply:	2 AAA batteries	
Display:	Body Fat Percentage:	4.0% to 50.0%
	Body Mass Index:	7.0 to 90.0
Set Ranges:	Height:	3' 4" to 6' 6" (100 cm to 199.5 cm)
	Weight:	23 lbs to 440 1/2 lbs (10 kg to 199.8 kg)
	Age:	NORMAL: 10 to 80 years old / ATHLETE: 18 to 60 years old
	Gender:	Male / Female
Operating Temperature and Humidity:	50°F to 104°F (+10°C to +40°C), 30% to 85% RH	
Storage Temperature and Humidity:	-4°F to +140°F (-20°C to +60°C), 10% to 95% RH	
External Dimensions:	Approx. 8" (L) x 5" (H) x 2" (W) (197 (L) x 128 (H) x 49 mm (W))	
Weight:	Approx. 8 oz (230 g) (not including batteries)	
Contents:	1 Body Fat Monitor, 2 AAA Batteries and 1 Instruction Manual	

Specifications and appearance are subject to change without notice for purposes of improvement.



OMRON HEALTHCARE, INC.

1200 Lakeside Drive
 Bannockburn, Illinois 60015
 Telephone: 847-680-6200
 Toll Free Number: 1-800-231-4030
 Fax Number: 1-800-637-6763
 Visit our Website at www.omronhealthcare.com

© 7/2009 Omron Healthcare, Inc.

Printed in the U.S.A.
 HBF-306CBLK PS Rev