

# **BASELINE<sup>®</sup>** **EVALUATION INSTRUMENTS**

## **12-inch 360° plastic goniometer 12-1000**

- permits observation of joint's axis of motion, and range of motion
- 360° head has three scales for ease of use
- primary scale reads in 1° increments
- goniometer arm has a linear scale in inches and centimeters
- goniometer arm remains at measurement until reset

---

### **other sizes and style plastic goniometers available**

12 inch 360° .....	12-1000*
8 inch 360° .....	12-1001*
6 inch 360° .....	12-1002*
6 inch 180° pocket.....	12-1005*
6 inch 360° rulongmeter.....	12-1006*
3 inch -30° to 180° finger.....	12-1012*

\*available standard or with white HiRes™ strip behind numbers (add HR to end of product code for HiRes™ model)



Fabrication Enterprises Inc.  
P.O. Box 1500  
White Plains, NY 10602 (USA)  
tel: 1-914-345-9300  
fax: 1-914-345-9800  
[www.Fab-Ent.com](http://www.Fab-Ent.com)



Authorized CE representative:  
RMS UK Ltd.  
28 Trinity Road  
Nailsea, Somerset BS48 4NU (UK)

